


60+ CLINICAL ARTICLES

More than 60 Thyroid RFA articles have been published with STARmed electrodes. These articles prove the success of this treatment by the reduction in thyroid nodule size, symptom and cosmetic scores, and Quality of Life.

* Product specification

Model name Specification	VIVA RF Electrode	star RF Electrode_Fixed	star Bipolar RF Electrode	star RF Electrode_Uni-Faced	star injectable RF Electrode
Gauge	18G	18/19G	17G	17G	18/19G
Length	70 -1 00 mm	70 -100 mm	70 - 100 mm	70 - 100 mm	70 - 100 mm
Active Tip	5 - 30 mm	3.8, 5, 7, 10, 15 mm	12, 14, 15, 16 mm	5, 7, 10, 15 mm	7, 10, 15 mm



www.STARmed4u.com
e-Mail : inquiry@STARmed4u.com
Tel : +82 506 816-3546 Fax : +82 506 816-4546
ST-CA-TH(I) Rev.0




Thyroid Radiofrequency ablation

Minimal Invasive Treatment for Thyroid Nodule



STARmed Electrodes

The sharpest electrode that needs no incision.
Clear visualization of electrode tip on ultrasound.
Easy & smooth electrode control for "Moving Shot" Technique.
Complete range of electrode type & specification.



Generator



Advanced Generator system
Dedicated algorithm to prevent charring.
Impedance sound guide for better control of procedure

About **STAR**med

Global Leader in Thyroid RFA

Largest number of published clinical studies done with STARmed devices
Partnership with Key Opinion Leaders for Thyroid RFA procedure training
New device and clinical procedure development based on customer feedback
Continuous professional educational activities held around the world

Globally Recognized Treatment Option

Pioneered Thyroid RFA into International Guidelines (Austria, Italy, Korea, U.K., U.S.A.)
Improves patient quality of life *

